

Consistency You Can Count On.

ELASTOMERICS

LEARN MORE



Consistency is key with Elastomerics.

From the beginning of a spool to the end, you need consistent force and reliable retention.



Made from medical-grade polymers

Strict quality control standards ensure material reliability.



Precision Manufacturing Process

Developed to minimize hydrophilic uptake and deliver consistent forces.



In-process testing and recalibration

Ensure every product meets the stringent force value tolerances you require.

LIGATURES

Made from **high-quality polymer** for optimal elasticity and engagement, our ligature ties are ergonomically designed for easy chairside application.

- Variety of sizes and feels
- Multiple colors and designs



Available in 33 colors

1 Bambino-Ties™

- Hold securely – or wear on fingers for ease of engagement and multi-color application
- 12 Ties per Stick, 84 Sticks/pkg (1,008 Ties total)

2 Versa-Ties™


- Score line allows easy conversion into half-sticks, helping reduce waste and cross-contamination.
- 22 Ties Per Stick, 46 Sticks/pkg (1,012 Ties total)

3 Select-Ties


- 20 ties per stick – eliminates waste and cross-contamination
- 20 Ties Per Stick, 50 Sticks/pkg (1,000 Ties total)

4 Mini-Ties



- Etched thumb tab – for secure grip with gloves
- 11 Ties Per Stick, 92 Sticks/pkg (1,012 Ties total)



TOP SELLER



TOP SELLER



INTRAORAL ELASTICS

Made from surgical-grade natural latex rubber, G&H intraoral elastics go through continuous testing, measuring and adjustments to ensure each batch **delivers the desired force every time**. Our elastics are available in a variety of colors and themes, as well as non-latex options.

Packs of 100 elastics—50 patient packs per box (5,000 total elastics).

LEARN MORE



SIZES	1/8" (3.2mm)	3/16" (4.8mm)	1/4" (6.4mm)	5/16" (7.9mm)	3/8" (9.5mm)
FORCES	LIGHT 2.5 oz.	LIGHT 2.5 oz.	LIGHT 2.5 oz.	LIGHT 2.5 oz.	LIGHT 2.5 oz.
	MEDIUM 4.5 oz.	MEDIUM 4.5 oz.	MEDIUM 4.5 oz.	MEDIUM 4.5 oz.	MEDIUM 4.5 oz.
	HEAVY 6.5 oz.	HEAVY 6.5 oz.	HEAVY 6.5 oz.	HEAVY 6.5 oz.	HEAVY 6.5 oz.
		EX. HEAVY 8.5 oz.	EX. HEAVY 8.5 oz.	EX. HEAVY 8.5 oz.	

Available in amber and assorted colors. Options vary by force and size.




FORCES	LIGHT 3.5 oz.	MEDIUM 4.5 oz.	HEAVY 6.0 oz.
SIZES		1/8" (3.2mm)	
		3/16" (4.8mm)	
		1/4" (6.4mm)	
		5/16" (7.9mm)	

Available in amber option.

NON-LATEX ELASTICS

Ideal for patients with latex sensitivity, it is manufactured from a **unique polymer** that **eliminates the proteins and other potentially allergenic compounds** associated with latex.

- Available Natural Color
- Forces: Medium and Heavy
- Sizes: 3/16" and 1/4"
- Each bag contains 100 patient packs of 100 elastics (1,000 elastics total)



POWER CHAINS

Precision-manufactured and scientifically tested. We control production every step of the way to deliver **optimal elasticity and retention**. From the first piece of chain to the last, our in-house quality checks ensure you get continuous, consistent force and predictable tooth movement.



Dyna-Link™

High elasticity with superior memory and rebound

- Latex-Free
- 15 ft. spool
- Thickness: .026"
- Available in 33 colors



PowerLinX™

Retains up to 85% of it's original pull strength even after 6 weeks of engagement


- 6.56 ft. spool
- Thickness: .020"
- Available in clear or grey

Short Long Continuous Short Long

POWER THREADS

Made from latex-free, medical-grade elastomers, G&H Elastic Power Threads are **precisely extruded** for use in **rotations and traction** applications.

- 25 ft. spool
- Thickness: .025"/.030"
- Solid, Solid Slip-Not, Hollow
- Available in clear or grey



SEPARATORS

Provide **quick, dependable and efficient separation**. Made from medical-grade polymer and featuring high strength and memory, they are available in standard and radio opaque options.

Radiopaque Separators

- Easily identified on x-ray
- Blue color
- Small (.155" O.D. and .042" thick)
- Large (.193" O.D. and .049" thick)
- 1,000 per pack

Molded Blue Separators

- Ideal for posterior separator
- .180" in diameter
- 960 per pack

Saf Separators

- Regular(Blue) and Heavy (Purple)
- 960 per pack

